# beurer

## BM 96 Cardio

<b>DE</b> Anhang für den behandelnden Arzt BM 96 EKG-Stick ..... 2	<b>IT</b> Appendice per il medico curante BM 96 stick ECG ..... 14
<b>EN</b> Annex for attending doctor BM 96 ECG stick ..... 5	<b>TR</b> Tedaviyi yürüten doktor için ek BM 96 EKG Çubuğu ..... 17
<b>FR</b> Annexe pour le médecin traitant Stick ECG BM 96 ..... 8	<b>RU</b> Приложение для лечащего врача Электрокардиограф BM 96 ..... 20
<b>ES</b> Apéndice para el médico Dispositivo de ECG BM 96 ..... 11	<b>PL</b> Załącznik dla lekarza prowadzącego BM 96 Element EKG ..... 23

## DEUTSCH

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## ENGLISH

The ECG stick from the BM 96 range is a single-channel ECG recording device for personal use. It records the ECG in time intervals of 30-120 [s]. The electrodes are integrated into the device. Through their positioning, they are oriented such that the results approximately show frontal leads according to Einthoven. The following applies in this context:

Method A: approx. lead II

Method B: approx. lead III

Method C: approx. lead I

Measured electrocardiograms are recorded using the ECG stick and can then be transferred to smartphones and computers on which the corresponding “beurer HealthManager Pro”/“beurer HealthManager App” applications have been installed. All measurements are available to the user and their attending doctor in both on-screen and printout format.

All displayed electrocardiograms are scaled by amplitude and time (recording speed: user can choose between 25 [mm/s] and 50 [mm/s] to display the results. Amplitude scaling: user can choose between 5 [mm/mV], 10 [mm/mV] and 20 [mm/mV] to display the results). The ECG printout enables the user to document basic patient data as well as information on the current ECG recording (e.g. current complaints/medication/the selected lead, etc.). This information can be passed on to the attending doctor to assist with the assessment of the ECG recording.


### **The described presentation of the ECG recording allows**

- (1) a doctor to subsequently provide an independent assessment of the results
- (2) for an automated, simplified assessment of the ECG results by the BM 96

The automated assessment mentioned under point (2) is not a complete or detailed clinical ECG diagnosis, but rather merely provides an initial starting point for the user due to the simplified nature of the data.

The status messages from the 30-120 [s] single-channel ECG recording are not a substitute for a medical diagnosis of any abnormal results.

**The BM 96 status messages are:**

	ECG recording is OK.
	Indicates a decreased heart rate (bradycardia) of less than 60 [bpm].
	Indicates an increased heart rate (tachycardia) of greater than 100 [bpm].
	Indicates atrial fibrillation

**The BM 96 device is not suitable for clinical emergency diagnoses.** Status messages do not inform of the presence of life-threatening disruptions to the heart rhythm such as ventricular fibrillation and/or asystole.

Status messages also do not show fundamental conduction disturbances of the heart. This applies to atrioventricular block patterns as well as complete or partial bundle branch block patterns on the left or right sides.

The primary function of the BM 96 device is to enable the recording of an electrocardiogram at home, with the aim of providing a simple means of documenting changes that could occur in the electrocardiogram in the absence of medical diagnostics. It is possible that home users may incorrectly record ECGs, e.g. due to motion artefacts, incorrect positioning of the electrodes or insufficient electrical contact (the presence of marginally smaller ECG amplitudes is signalled). As a result, any subsequent ECG diagnoses may be more difficult or even prevented.

We recommend that users undergo a medical examination if pathological status messages repeatedly occur, in particular.

**Regardless of the results delivered by the BM 96, it is always necessary for users to immediately undergo a medical examination if they experience cardiac symptoms that could indicate acute medical problems.**

Whilst the recording of BM 96 electrocardiograms can never replace a medical cardiological diagnosis, the submission of electrocardiogram results may assist with the diagnosis.

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BEURER GmbH • Söflinger Str. 218 • 89077 Ulm (Germany)